

T AND EXPORT

Anesthesia and Respiratory

Aerosol Chamber 储雾罐

Function

- Aerosol Chamber is used to help deliver medication to cats/dogs suffering from respiratory conditions such as chronic bronchitis, laryngeal paralysis or tracheal collapse.

Material

- Antistatic Chamber is PP.
- Silicone connector.
- Soft silicone mask, 3 sizes available.
- Safety silicone breathing valve.

Certificate

- ISO/SGS/Antistatic report.



Aerosol Chamber Item No.	Mask Item No.	Mask Size	Mask OD	Fit Weight
SAC-5	C3-0	0#	51.1mm	0-5kg
	C3-1	1#	63.9mm	5-10kg
	C3-2	2#	78.5mm	>10kg

Color can be customized

- Green, Yellow, Pink 3 colors available.



Taking care of animal health...