

TANDEXPORT





Cohesive Bandage(Latex / Latex Free)

Size: 2.5cm*4.5m / 5cm*4.5m / 7.5cm*4.5m / 10cm*4.5m / 15cm*4.5m Non-standard specifications can be customized.

Color: All sorts of color, Can print a variety of patterns and LOGO.

Cloth: Non-woven fabric.

Glue: Latex / Latex Free.

Purpose: 1. Medical fixation, orthopedic splint fixation, wound dressing fixation, catheter fixation;

2. Joint protection before exercise;

3. Protect students' fingers, nail protection, tattoo protection;

4. Animal medical fixation, animal decoration.

ANIMAL



01 奶牛
Cow



02 长颈鹿
Giraffe



03 豹纹
Leopard



04 蟒蛇
Snake



05 虎纹
Tiger



01 斑马纹
Zebra

CAMO



01 丛林迷彩
Jungle Camo



02 军绿迷彩
Olive Camo



03 树皮迷彩
Wood Camo



04 雪地迷彩
Snow Camo



05 沙漠迷彩
Desert Camo



06 ACU 迷彩
ACU Camo

DARK



01 骷髅
Skull



02 爱
Love



03 洛丽塔
Lolita



04 闪电
Lightning

COLORFUL



01 樱桃
Cherry



02 柠檬
Lemon



03 草莓
Strawberry



04 蓝莓
Blueberry



05 牛油果
Avocado



06 水蜜桃
Peaches



07 黄泥笑脸
Yellow smile



08 白底爱心
Pink Heart



09 熊猫
Panda



Cotton flexible cohesive bandage (Latex)

Size: 2.5cm*4.5m / 5cm*4.5m / 7.5cm*4.5m / 10cm*4.5m / 15cm*4.5m Non-standard specifications can be customized.

Color: All sorts of color, Can print a variety of patterns and LOGO.

Cloth: Cotton.

Glue: Latex.

Purpose: 1. Medical fixation, orthopedic splint fixation, wound dressing fixation, catheter fixation;
2. Joint protection before exercise.



Cohesive gauze

Size: 2.5cm*4.5m / 2.5cm*9m / 1.8cm*18m Non-standard specifications can be customized.

Color: All sorts of color.

Cloth: Gauze.

Glue: Latex / Latex Free.

Purpose: 1. Medical fixation, orthopedic splint fixation, wound dressing fixation, catheter fixation;
2. Protect fingers.



PBT Cohesive bandage

Size: 2.5cm*4.5m / 6cm*4.5m / 8cm*4.5m / 10cm*4.5m Non-standard specifications can be customized.

Color: All sorts of color.

Cloth: PBT Gauze.

Glue: Latex/Latex Free.

Purpose: Medical fixation, orthopedic splint fixation, wound dressing fixation, catheter fixation;



Foam bandage

Size: 6cm*1m / 6cm*5m.

Color: All sorts of color.

Cloth: Foam.

Glue: Latex /Latex Free.

Purpose: Protect fingers.



Flat edge sports tape / Serrated sports tape
Perforated athletic tape
LockerTape rayon micro porous rigid sports tape

Size: 2.5cm*13.7m / 3.8cm*13.7m / 5cm*13.7m / 2.5cm*9.1m / 3.8cm*9.1m / 5cm*9.1m Non-standard specifications can be customized.

Color: White and Color, can print a variety of patterns and LOGO.

Cloth: Cotton/terylene/rayon.

Glue: Hot melt /zinc oxide.

Purpose: Used for pre-exercise protection.





HockeyTape

Size: 2.5cm*10m / 2.5cm*20m / 2.5cm*25m / 3.8cm*20m / 3.8cm*25m Non-standard specifications can be customized.

Color: White and Color, can print a variety of patterns and LOGO.

Glue: Hot melt.

Purpose: Used to attach hockey sticks.



Underwrap

Size: 7cm*20m / 7cm*27m / 6.7cm*27.4m Non-standard specifications can be customized.

Color: White and Color, can print a variety of patterns and LOGO.

Cloth: Foam.

Purpose: protect the skin of athletes.



PE Tape

Size: 2.5cm*9.1m / 5cm*9.1m / 7.5cm*9.1m Non-standard specifications can be customized.

Color: Transparent.

Glue: Hot melt / acrylic adhesive.

Purpose: 1. Medical pipeline fixation, infusion fixation.
2. Beauty eyelash.



Non-woven paper tape

Size: 2.5cm*9.1m / 5cm*9.1m / 7.5cm*9.1m Non-standard specifications can be customized.

Color: White and Skin.

Glue: Hot melt.

Purpose: Used for medical pipeline fixation, infusion fixation.




Non-woven tape

Size: 2.5cm*9.1m / 5cm*9.1m / 7.5cm*9.1m / 10cm*9.1m Non-standard specifications can be customized.

Color: White and Skin.

Glue: Hot melt /acrylic adhesive.

Purpose:Used for medical pipeline fixation, infusion fixation.



Silk tape

Size: 2.5cm*9.1m / 5cm*9.1m / 7.5cm*9.1m / 10cm*9.1m Non-standard specifications can be customized.

Color: White.

Glue: Hot melt /acrylic adhesive.

Purpose:Used for medical pipeline fixation, infusion fixation.



Plastic Zinc oxide tape

Size: 2.5cm*9.1m / 5cm*9.1m / 7.5cm*9.1m / 10cm*9.1m Non-standard specifications can be customized.

Color: White and Skin.

Glue: Zinc oxide hot melt adhesive.

Purpose: Used for medical pipeline fixation, infusion fixation.




Silk tape

Size: 1.25cm*5m / 2.5cm*5m / 5cm*5m / 7.5cm*5m Non-standard specifications can be customized.

Color: White and Skin.

Glue: Zinc oxide hot melt adhesive.

Purpose: Used for medical pipeline fixation, infusion fixation.



Foam tape

Size: 2.5cm*4.5m / 2.5cm*5m.

Color: Skin.

Glue: Acrylic adhesive.

Purpose: 1. Used for medical fixation For heel anti-wear care;
2. For heel anti-wear care.




Foam tape

Size: 2.5cm*4.5m / 2.5cm*5m.

Color: White.

Glue: Acrylic adhesive.

Purpose: 1. Used for medical fixation For heel anti-wear care;
2. For heel anti-wear care.




Easy tear tape

Size: 1.25cm*10m / 2.5cm*10m / 1.25cm*5m / 2.5cm*5m Non-standard specifications can be customized.

Color: A variety of color

Glue: Hot melt /acrylic adhesive.

Purpose:1.Used for medical pipeline fixation, infusion fixation ;
2.Beauty eyelash.



Non-woven dressing roll

Size: 宽度(2.5cm-10cm) * (5-20)m Non-standard specifications can be customized.

Color: White.

Glue: Hot melt /acrylic adhesive.

Purpose:Used for medical pipeline fixation, infusion fixation, dressing fixation.



Kinesiology tape

Size: 2.5cm*5m / 3.8cm*5m / 5cm*5m / 5cm*31m Non-standard specifications can be customized.

Color: All sorts of color, Can print a variety of patterns and LOGO.

Glue: Hot melt / acrylic adhesive.

Purpose: 1. Sports protection, joint protection, prevent muscle strain, improve exercise function

2. Physical therapy.





Die-cut kinesiology tape

Size: 2.5cm*5m / 3.8cm*5m / 5cm*5m / 5cm*31m Non-standard specifications can be customized.

Color: All sorts of color, Can print a variety of patterns and LOGO.

Glue: Hot melt /acrylic adhesive.

Purpose: 1.Sports protection, joint protection, prevent muscle strain, improve exercise function
2.Physical therapy.



Shape pre-cut kinesiology tape

Size: 2.5cm*10cm / 5cm*10cm / 5cm*25cm Non-standard specifications can be customized.

Color: All sorts of color, Can print a variety of patterns and LOGO.

Glue: Hot melt /acrylic adhesive.

Purpose: 1.Sports protection, joint protection, prevent muscle strain, improve exercise function
2.Physical therapy.




Boob Tape

Size: 5cm*5m.

Color: All sorts of color, Can print a variety of patterns and LOGO.

Glue: Hot melt.

Purpose: Used for female breast lifting and shaping.




Soft close skin Boob Tape

Size: 5cm*5m Non-standardspecifications can be customized.

Color: All sorts of color, Can print a variety of patterns and LOGO.

Glue: silica gel.

Purpose: It is used for female breast lifting and shaping; It can be used repeatedly.



Cross tape

Size: A2 patch / A4 patch / A6 patch / A9 patch, Non-standard specifications can be customized.

Color: All sorts of color.

Glue: Acrylic adhesive.

Purpose: 1. Used for acupuncture and physiotherapy fixation;

2. It is used to prevent swelling and itching after mosquito bites.



Elastic adhesive bandage

Size: 2.5cm*4.5m / 5cm*4.5m / 7.5cm*4.5m / 10cm*4.5m Non-standard specifications can be customized.

Color: White/Skin.

Glue: Hot melt / zinc oxide.

Purpose: 1. Used for medical fixation;

2. Heavy duty exercise protection and joint fixation.



LightWrap Tear-Light elastic adhesive bandage

Size: 2.5cm*4.5m / 5cm*4.5m / 7.5cm*4.5m / 10cm*4.5m Non-standardspecifications can be customized..

Color: All sorts of color, Can print a variety of patterns and LOGO.

Glue: Hot melt.

Purpose:1.Used for medical fixation;

2.Light sports protection and joint fixation.



Skin tone high elastic bandage

Size: 2.5cm*4.5m / 5cm*4.5m / 7.5cm*4.5m / 10cm*4.5m Non-standardspecifications can be customized.

Color: All sorts of color.

Glue: Hot melt.

Purpose:1.Used for medical fixation;

2.Heavy duty exercise protection and joint fixation.



LightWrap Tear-Light elastic adhesive bandage

Size: 2.5cm*4.5m / 5cm*4.5m / 7.5cm*4.5m / 10cm*4.5m Non-standardspecifications can be customized.

Color: All sorts of color, Can print avariety of pattemns and LOGO.

Glue: Hot melt.

Purpose:1.Used for medical fixation;

2.Heavy duty exercise protection and joint fixation.

