






## Cat Catheter 猫用导尿管

- Made of medical grade nylon materia

| Item No. | Size(OD) | Length(mm) | Without stylet/Hole types   |
|----------|----------|------------|---|
| PCC-1    | 1.0mm    | 130        |  end open  |
| PCC-2    | 1.3mm    | 130        |  side open |





| Item No. | Size(OD)       | Length(mm) | With stylet/Hole types  |
|----------|----------------|------------|---|
| PCCT-1   | 3Fr endopen    | 120        |  end open    |
| PCCT-2   | 3Fr closed end | 120        |  closed end |



## Silicone Cat Catheter 硅胶猫用导尿管

### Function

- Made of imported medical grade silicone.
- Soft tube is more comfortable and can relieve the pain of the patient.
- Remark: size and length can be customized.

| Item No. | Size(Fr) | Length(mm) | With stylet/Hole types  |
|----------|----------|------------|---|
| SCCT-1   | 3        | 130        |  end open  |
|          |          |            |  side open |



Taking care of animal health...