

T AND EXPORT

Anesthesia and Respiratory

Endotracheal Tube 气管插管

- Made of medical grade PVC, clear and smooth.
- High volume, low pressure cuff maintains good sealing.
- Murphy eye design, tube with X-ray opaque line.

Notes

- Uncuffed type is available, from ID2.0-ID9.5.



Item No.	Size ID (mm)	Item No.	Size ID (mm)
ET25PC	2.5	ET75PC	7.5
ET30PC	3.0	ET80PC	8.0
ET35PC	3.5	ET85PC	8.5
ET40PC	4.0	ET90PC	9.0
ET45PC	4.5	ET95PC	9.5
ET50PC	5.0	ET100PC	10.0
ET55PC	5.5	ET105PC	10.5
ET60PC	6.0	ET110PC	11.0
ET65PC	6.5	ET120PC	12.0
ET70PC	7.0	ET130PC	13.0

Reinforced Endotracheal Tube 加强气管插管

- Made of medical grade PVC, clear and smooth.
- High volume, low pressure cuff maintains good sealing.
- Tube with spiral reinforced wire reducing risk of kinking.
- Murphy eye design, tube with X-ray opaque line.

Notes

- Can be equipped with stylet.
- Uncuffed type is from ID2.0-ID9.5



Item No.	Size ID (mm)	Item No.	Size ID (mm)
RET25PC	2.5	RET65PC	6.5
RET30PC	3.0	RET70PC	7.0
RET35PC	3.5	RET75PC	7.5
RET40PC	4.0	RET80PC	8.0
RET45PC	4.5	RET85PC	8.5
RET50PC	5.0	RET90PC	9.0
RET55PC	5.5	RET95PC	9.5
RET60PC	6.0	RET100PC	10.0

Taking care of animal health...