
JX7900Aplus Large Animal Anesthesia Machine (JX7900A+JX560G)



JX7900Aplus model is an excellent large animal anesthesia machine. It is specially used for large animal, like horse, cow, tiger, pig etc. It can be used in farm, zoo, wild animal zoo, forest park, animal research center and so on.

Large Animal Anesthesia Machine(JX7900A)

Technology Specification:

Inspiration mode: open mode, closed mode, semi-open and semi-closed mode

Gas source: O₂, 0.3--0.4Mpa

O₂ Flush valve: 25L/min--75L/min

Flow meter: high precision O₂ flow meter, 0 --10L/min, O₂&Air two tubes flow meter is option.

CO₂ circle absorber: with soda lime chamber, Pop off valve, airway pressure gauge, inspiration valve and expiration valve, the inspiration and expiration valve is demo-type and visible.

Airway pressure gauge: -20cmH₂O--100cmH₂O &-2Kpa--10Kpa

Soda lime chamber: 8Liter, it is easy to change soda lime and used for over 1000kg animal.

Bains connector: it is used for open breath circuit, connect with Jackson breath circuit or baird breath circuit.

Anesthesia vaporizer: TEC3 vaporizer. Enflurane, Isoflurane, Halothane, Sevoflurane, one vaporizer with selectatec compatible mounting bar.

A handle on the top

Large Animal Anesthesia Ventilator(JX560G)

Technology specification:

Driving type: Gas Driving Electronic Control

Driving gas: O₂ or Air, 0.35Mpa-0.6Mpa, max flow more than 400l/Min at 0.35Mpa

Breath Mode: IPPV

Inspiration time: 0.2-10S

Breath Rate: 2-80

Tidal Volume: 0-20L

Airway Pressure Gauge Monitor: 0-99cmH₂O

PEEP: when PEEP is more than 5cmH₂O, Pop Off valve will be open

Bellow: 20L

Inner Battery: it will work 3 hours

Power: AC100-240V, 50-60HZ

Alarm Parameters:

Gas source alarm: when gas source pressure is less than 0.35Map, it will alarm;

AC power alarm: when AC power is faulty, it will alarm;

IRV alarm: when I:E is IRV ventilation, it will alarm;

Low Battery: when inner batter is low, it will alarm;

Standard configuration:

Anesthesia machine, anesthesia ventilator, one unit X60 anesthesia vaporizer, one set breath circuit, one unit 15L airbag, one unit O2 gas source pipeline, one unit Power cable, trolley.

Packing information:

Including Two paper cartons:

No1 carton: JX7900A Large Animal Anesthesia Machine with X60 anesthesia vaporizer:62*62*62cm, gross weight:46kg

No2.carton: JX560G Large Animal Anesthesia Ventilator(including trolley): 65cm*64cm*143cm, gross weight:50kg.