


生产实施 »

» PRODUCTION EQUIPMENT



检验实施 »

» INSPECTION EQUIPMENT



产品介绍 »

» PRODUCT INTRODUCTION



Pre-powdered Sterile Latex Surgical Gloves



PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Lightly powdered, powdered content≤10mg/dm ²
Color	Natural
Thickness	Standard
Cuff Style	Beaded
Mould Style	Anatomically curved fingers
External Glove Surface	Palm Textured & other place Smooth
Internal Glove Surface	Polymer Coated
Sterile	Yes, Sterilized by Rays
Freedom from Holes (Inspection level G-I)	1.0 AQL
Product Standards	ASTM D3577, EN 455 1-4, ISO 10282
Quality Standards	ISO 13485:2016
Packaging	50 pairs per box, 10 boxes per carton, 1040 cartons per 40HQ container 450 cartons per 20GP container 50 pairs per box, 8 boxes per carton, 1300 cartons per 40HQ container 550 cartons per 20GP container
Shelf Life	5 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

PRODUCTS CODES INFORMATION

Size	6	6.5	7	7.5	8	8.5
Product Code	9101	9102	9103	9104	9105	9106

产品介绍 »

» PRODUCT INTRODUCTION



Powder free Sterile Latex Surgical Gloves



PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder free
Color	Natural
Thickness	Standard
Cuff Style	Beaded
Mould Style	Anatomically curved fingers
External Glove Surface	Palm Textured & other place Smooth
Internal Glove Surface	Polymer Coated
Sterile	Yes, Sterilized by Rays
Freedom from Holes (Inspection level G-I)	1.0 AQL
Product Standards	ASTM D3577, EN 455 1-4, ISO 10282
Quality Standards	ISO 13485:2016
Packaging	50 pairs per box, 10 boxes per carton, 1040 cartons per 40HQ container 450 cartons per 20GP container 50 pairs per box, 8 boxes per carton, 1300 cartons per 40HQ container 550 cartons per 20GP container
Shelf Life	5 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

PRODUCTS CODES INFORMATION

Size	6	6.5	7	7.5	8	8.5
Product Code	9108	9109	9110	9111	9112	9113

产品介绍 »

» PRODUCT INTRODUCTION



Powder free Sterile
Latex Surgical Gloves



PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder free
Color	Natural
Thickness	Thickened (Smooth area: Min 0.14mm, Textured area: Min 0.2mm)
Cuff Style	Beaded
Mould Style	Anatomically curved fingers
External Glove Surface	Palm Textured & other place Smooth
Internal Glove Surface	Polymer Coated
Sterile	Yes, Sterilized by Rays
Freedom from Holes (Inspection level G-I)	0.65 AQL
Product Standards	ASTM D3577, EN 455 1-4, ISO 10282
Quality Standards	ISO 13485:2016
Packaging	50 pairs per box, 10 boxes per carton, 960 cartons per 40HQ container 390 cartons per 20GP container 50 pairs per box, 8 boxes per carton, 1180 cartons per 40HQ container 480 cartons per 20GP container
Shelf Life	5 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

PRODUCTS CODES INFORMATION

Size	6	6.5	7	7.5	8	8.5
Product Code	9115	9116	9117	9118	9119	9120