

T AND EXPORT

ELECTRICAL STIMULATION EJACULATOR

The portable electrical stimulation ejaculator can achieve the collection of animal semen by electrically stimulating animals and inducing their ejaculation behavior. This product is easy to operate, safe and reliable, and has high collection efficiency. It is widely used in animal husbandry, especially for forced semen collection from animals that cannot be domesticated, have no sexual desire or are too wild males.

Neutral products are available and can be customized (logo, packaging, instruction manual).



Usage scenario



Product parameters

Stimulation voltage: 2.5 - 25V
Stimulation frequency: 5Hz - 100Hz
Automatic stimulation rhythm:
0.5 times/second - 2.5 times/second
Charging time: 3 - 4 hours

Stimulating rods for cattle:
Probe length: 411mm, diameter: 25mm
Stimulating rods for sheep:
Probe length: 232.5mm, diameter: 20mm
Stimulating rods for rabbits:
Probe length: 115mm, diameter: 6mm
Stimulating rods for mice:
Probe length: 40mm, diameter: 2mm

Packaging dimensions



Dimensions: 30*21*16cm
Weight: 2.3kg
(with sheep-shaped tickler)

Probe Introductions



Tortoise, Mouse



Monkey, Fox,
Cat, Canine



Pig, Sheep/Goat, Deer



Cattle, Horse,
Donkey, Camel

Product display



Recommended combination



Pig Semen Collection
Throstatic Cup



Semen Collection Thawed
Throstatic Cup



Sperm Quality Analyzer



Convenient Microscope



constant temperature
microscope



Pig Pregnancy Test Kit