

S150 Portable Slit Lamp



Feature

- LED illumination:
 - Magnification: 5x
 - Comfortable using experience:
 - Power saving: S150 can work for long hours without changing batteries.
- Can be used for humans and animals



Specifications:

Magnification:	5x
Length of slit:	5mm ~ 11mm (stepless)
Minimum width of slit:	0.2mm
Working distance:	12 ~ 32mm
Light source:	LED (3w/1w, 20,000 hours working life)
Power supply:	AA alkaline batteries x 2 (more than 3 hrs working life)
Net weight:	150g(without batteries)
Guarantee:	1 Year

YZ-9 Exophthalmometer



90D 78D 20D Sphere Lens



90D



78D



20D