

Veterinary Anesthesia Mask 动物面罩



- Face mask for dogs, cats and other small animals devised for anesthesia, ventilation, aerosol therapy.
- According to the animal- patient's face size, choose correct mask size.
- Made of a high transparency PC cone and a soft black silicone diaphragm.
- Available in 6 sizes, 1#, 2#, 3#, 4#, 5#, 6#.
- Both PC ring and stainless ring available
- Standard connectors (15mm, 22mm).



Warning

- One mask for one animal patient use only.

Item No.	Size	Hight	ID	OD	Connector
VM-1	1#	44mm	24mm	48mm	15mm
VM-2	2#	49mm	28mm	53mm	
VM-3	3#	58mm	35mm	64mm	22mm
VM-4	4#	98mm	29mm	92mm	
VM-5	5#	132mm	49mm	123mm	
VM-6	6#	141mm	71mm	136mm	



● Silicone Neck Strap



● PP Ring



● Stainless Ring

Taking care of animal health...